

## User interfaces


### ILL

See separate file


### ANSTO, Quokka SANS

NIST macros

<https://www.ansto.gov.au/sites/default/files/2022-03/QUOKKA%20user%20manual%20March%202022.pdf>



Leaving 'ask' will mean that during the reduction Igor will prompt you to select the file you wish to use:




ANSTO Bilby SANS, no GUI, files & reduction are kept in the csv files:

input_9862.csv															
index	T_EmptyBeam	T_Sample	Sample	thickness [cm]	mask_transmission	mask	BlockedBeam	suffix	additional_description	T_BlockedBeam	mask_transmission	StartTime	EndTime		
1															
2															
3															
4															
5	0	BBY0047802	BBY0047802	BBY0047804	0.2	mask_transmission_mask_scattering_9862_rearpackoff	BBY0047806	EmptyBeam							
6	1	BBY0047803	BBY0047803	BBY0047804	0.2	mask_transmission_mask_scattering_9862_rearpackoff	BBY0047806	EmptyBeam							

settings_9862.csv															
index	csv_file_name	reduced_files_folder	binning_wavelength_in	binning_q	binning_wavelength_trai	RadiusCut	WaveCut	transmission_fit	PolynomialOrder	wavelength_intervals	wav_delta	reduce_2D	2D_number_data_points	plot_2D	
1															
2															
3															
4															
5															
6	0	input_9862.csv	data	2.0, -0.1, 20.0	0.0005, -0.05, 1.5	2.0, -0.1, 20.0	0	0	Polynomial	3	TRUE	2	FALSE	100	FALSE
7															
8	1	input_9862.csv	data_final	3.0, -0.1, 16.0	0.0005, -0.1, 1.5	2.0, -0.1, 20.0	0	0	Polynomial	3	FALSE	2	FALSE	100	FALSE
9															

**ANSTO ToF spectrometer Pelican** (this is the current state, shall be improved) (to note: backscattering Emu has similar User interface; Emu reduction code is still not in the release):



There is a QENS fitting tool that ISIS developed for Mantid that seems very useful that we are only just starting to use.

I also want to learn about the VATES option for single crystal processing (as currently we split this into two and start with Mantid, but finish with Horace/matlab).

## CSNS

The attached figure is the batch processing interface we use currently. RunNo for each sample and its transmission and EB/EC need to be type in manually, and by clicking add button, these information will be stored in the spreadsheet on the bottom, and can be edited in the spreadsheet later. By clicking the CalAll button, every entry in the spreadsheet will be calculated sequentially.

We are hoping that the Batch processing interface could allow the upload of an excel or cvs file, from where the relevant run number and other information for each sample can be extracted directly.

SANS Data Reduction V1.2.19

Configuration Transmission Scattering IQChecking **Batch** IQ2D CalibrationRatio Anisotropic Slice

Transmission template

Empty Trans Run: 14119

Trans Fit Type: Raw

Wavelength(AA): Min: Step: Max:

1 0.01 9.8

Scattering template

Empty Direct Run: 14409

Q: Min: Step: Max:

0.051 -0.08 1.2 AutoQRange

Sample info

Sample Trans Run: 14200

Sample Scatt Run: 14219

Thickness: 1.14 mm

Scaling Factor: 1.0

Add Remove RemoveSampleNo: 3

	EmptyTrans	SampleTrans	EmptycellTrans	Emptycell-EmptyTrans	SolventTrans	Solvent-EmptyTrans
1	14199	14200	NA	NA	NA	NA
2	14199	14200	NA	NA	NA	NA
3	14199	14200	NA	NA	NA	NA
4	14199	14200	NA	NA	NA	NA
5	14199	14200	NA	NA	NA	NA

CalAll Back Next

## GISANS

### ISIS

The GISANS is not in the priority list of the Large Scale Structure group so won't get done for a while so anything that has been done already and could be expanded in a collaborative way would be great.

### CSNS

In terms of GISANS, we hope Mantid could provide a specified GISANS data visualization interface similar to the current Instrument interface. We hope the interface would allow us to bin the data within certain wavelength range and make line cuts in vertical and horizontal directions and obtain and store the intensity distribution along the line (like imageJ).

The transmission correction for GISANS data has always been a problem for us as spallation neutron source. I am not sure if I could count on Mantid to solve this issue though. I would be grateful if we could discuss this issue in the further discussion

### ANSTO

I've talk to our scientist from our second SANS machine, they have GISANS set-up. For the data, they've used Grasp, just to get an integrated intensity across a chosen area of detectors. And I cannot see a wish to put an effort into the Mantid GISANS part.

## Algorithms that need improvements and validation tests

- Pascal Manuel:

Absorption corrections

- Anna Sokolova:

wide\_angle\_correction

(<https://docs.mantidproject.org/nightly/algorithms/SANSWideAngleCorrection-v1.html>).